



ORDER OF STEPS FOR NUTSCHE FILTER:

1. PRECIPITATION SUSPENSION (STREAM 205) ENTERS
2. NITROGEN ENTERS, AND THE TANK IS PRESSURIZED TO 200kPa DURING FILTERING
3. FILTRATE LEAVES TO TANK T-208
4. WASH SOLVENT (STREAM 207) ENTERS
5. WASH SOLVENT LEAVES TO TANK T-208
6. DISSOLUTION SOLVENT (STREAM 212) ENTERS
7. THE DISSOLVED SOLUTION (STREAM 208) EXITS TO STAGE 3





Stream Identification No.	309	310									NOTES:  CWS - Cool water supply CWR - Cool water return  SS - Steam supply CR - Condensate Return							
Stream Description																		
	Wt%	Wt%																
Component	Formula										LEGEND			REVISIONS				
Acetic Acid		0	0								<div>◇</div> Stream Number <div>○</div> Mass Flow (kg/hr) <div>□</div> Temperature (C) <div>○</div> Pressure (atm)	REV	DATE	DESCRIPTION				
Water		27.4	29.9							A		2/16/2022	First Draft of Process					
Acetonitrile		0	0							B		2/18/2022	Pre-Design Review					
Acetone		5.1	0							C		3/25/2022	Process Update					
Formic Acid		9.2	0							D		4/1/2022	Final Process Design					
Isopropanol		52.4	0															
Ciprofloxacin HCl		0.0173	70.1															
Tetrabutylammonium hydroxide		0	0															
Hydrochloric Acid		5.950	0															
trifluorobenzoyl)prop-2-enoate		0.006	0															
trifluorobenzoyl)acrylsaeure-ethylester		0.005	0															
1-cyclopropyl-6-fluoro-4-oxo-7-(piperazin-1-		0.006	0															
Total Mass Flow Rate (kg/batch)		57206194.3	23888.65								University of Michigan DEPARTMENT OF CHEMICAL ENGINEERING			COURSE NO. CHE 487	DRAWN BY ELIOT	DATE 3/25/2022		
Total Molar Flow Rate		1520.809476	0.605472											Group: 05		PROCESS FLOW DIAGRAM		
Max Pressure (atm)																		
Max Temperature											DRAWING TITLE			Crystallization - Stage 3			REV.	C

